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EN Intended use
For in vitro diagnostic use. The Control Cuvette Set is intended to be used for target value determination of a new Control Cuvette for the Hemo Control and the Hemo Vet photometer.

Product description
The Control Cuvette Set consists of a Control Cuvette, 2 Cleaner and the control solution Hb-con norm. The control solution contains 1ml ready-to-use solution, made of bovine hemolysate. The hemoglobin target values are determined by analyzes with the calibrated Hemo Control measuring system. The Cleaner is intended for single use only.

Important Notes for handling the Control Cuvette
Store the control cuvette always in the original box. Avoid contamination of the control cuvette with potential infectious material. Clean and disinfect the cuvette holder and the instrument before using the control cuvette. Keep away humidity (air humidity >85%) from the control cuvette. Do not clean the control cuvette with any solutions. If necessary clean the control cuvette carefully with a dry cotton bud.

Instructions for use
Clean the optical unit of the device and the cuvette holder.

- Remove the cuvette holder according to the user manual of the device.
- Take a Cleaner out of the packaging. Carefully insert the Cleaner into the slot of the cuvette holder until you feel a slight resistance.
- Slightly push down the Cleaner and move it further into the slot. Do not use any force, only slight resistance should be felt.
- Wipe frequently for cleaning the optical unit. Push the Cleaner up and down while wiping.
- Remove the Cleaner. Repeat the procedure with a new Cleaner if the used Cleaner is soiled.

- Clean the cuvette holder according to the user manual of the device.
- Wait at least 15 minutes before you reinsert the dry cuvette holder.

Check the measurement accuracy of the device.

- Check if an Adaptation is configured in the device (device menu - info, refer to the user manual of the device). The set parameters have to be considered (see point g).
- Take out the control solution from the fridge and allow it to warm up to room temperature for at least 15 minutes. Mix the control carefully by repeated rotating and rolling. Note the opening date on the bottle when opening the first time.
- Discard the first drop to avoid measurement errors due to incomplete mixture. Dispense a sufficiently large drop of solution (about 15µL) on a non-absorbent material (e.g. PE film).
- Fill and measure the Hemoglobin Microcuvette according to the user manual of the photometer. **Discard the first drop from the bottle at every use. Never do fill the cuvette directly from the bottle. Clean the bottle and the cap before closing.**
- Repeat the steps b) to d) five times in all and note the measurement results.
- Clean the cuvette holder and the device according to the user manual.
- IF ADAPTATION IS CONFIGURED: Calculate the real value x from the measured value y using the following equation: $x = (y - n) / m$
x...real value, y...measured value, n...offset, m...adaption factor

The expected results are given on the packaging. The measured results must fall within the specified limits.

If one result is outside the limits, check the correct execution of all procedure steps and the compliance with storage conditions and expiry dates of the control solution and the Microcuvettes. In case of doubt repeat measurement with a new bottle of control solution or Microcuvettes taken from a new can. If the results are outside the limits again, the device has to be checked by authorized service.

- If the results are within the limits, the control cuvette can be calibrated.
- Note the serial number of the device on the box of the control cuvette. Each control cuvette with indicated limit values belongs to a specific instrument.
 - Measure the control cuvette five times on the photometer and note the results.
 - Calculate the mean of the five results.
 - Round the value for mmol/l and g/dl units with one position after decimal point and for g/l unit as integer.
 - Note the determined value as target value on the box of the control cuvette in the designated lines.

Warning
Control solutions should be handled as potential infectious material, wear suitable gloves.

Disposal
The cleaner and the control solution must be disposed of after use as potentially infectious waste according to the regulations applicable to your establishment. Control cuvettes can be disposed of as residual waste.

Details		
	Reference number	3000-6138 Control Cuvette Set
	Lot number	see imprint
	Storage temperature	Hb-con norm: 2°C - 8°C Cleaner: 2°C - 30°C Control cuvette: 0°C - 40°C
	Expiry date (unopened)	see imprint
	Expiry date after opening	Hb-con norm: 30 days at 2°C - 30°C Cleaner: 10 min
	Opening date	note date

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